



Home (/) >> PRODUCT (<http://www.hssdgroup.com/PRODUCT/>) >> Metallographic Equipment (http://www.hssdgroup.com/PRODUCT/Metallographic_Equipment/)



MP-1C Metallographic Specimen Grinding & Polishing Machine

Product description:

1. Application and Feature: MP-1C Metallographic Specimen Grinding Polishing Machine is equipped with single disc and features stepless speed changing grinding and polishing. The rotating speed of grin



(http://www.hssdgroup.com/PRODUCT/Universal_Testing_Machine) (http://www.hssdgroup.com/PRODUCT/Universal_Testing_Machine) (http://www.hssdgroup.com/PRODUCT/Universal_Testing_Machine) (http://www.hssdgroup.com/PRODUCT/Universal_Testing_Machine) (http://www.hssdgroup.com/PRODUCT/Universal_Testing_Machine)

Share With:

We're here to help:
Easy ways to get the answers you need.

Send E-mail (<mailto:Louis@hssdtest.com>)

Chat Now (<skype:sales02@hssdtest.com?chat>)



Feature:

(/)

MP-1C Metallographic Specimen Grinding & Polishing Machine is equipped with single disc and features stepless speed changing grinding and polishing. The rotating speed of grinding and polishing disc is adjustable from 50 to 1000 rpm. It can realize the whole process of rough grinding, fine grinding, rough polishing and finishing polishing for specimen preparing. This machine features easy operation and economical applying. It is the ideal specimen preparing instrument for using in factories, scientific research institutes and laboratories of colleges.

2. Main Specifications:

Model	MP-1C
Grinding/Polishing Disc Diameter	200mm
Grinding Disc Rotating Speed	50-1000 rpm (stepless speed)
Motor	YSS6324 250 W
Dimension	400×730×310 mm
Weight	13 Kg
Operating Voltage	AC 220V,50Hz